

EasyMP™ beads (40-60 µm) in ethanol (95 vol%.) solution - different polymer types and colors available - please contact us for details

Art. ID MPS-013
Unit 30 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1170 III E1

Description

available colors: Clear/Red/Blue/Fluorescent green (Color only for PE)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
available polymer types	PE (polyethylene) / PMMA (Polymethyl methacrylate)		<1	%		