

EasyMP™ fibers (50-1000 µm) in ethanol (95 vol%.) solution - different polymer types and colors available - please contact us for details

Art. ID MPS-018
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1170 III E1

Description

available colors: Clear/Black/Red/Blue/Green/Violet

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
available polymer types	Polyester (PET) / Nylon (PA6,6) / Acrylic (PAN) / Rayon / Cotton / Si lk		<1	%		