

Custom Solution - 5 components in HNO3 5%

Art. ID LM24-CUS-73032
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Germanium (Ge)	[7440-56-4]	200	mg/L		
	Yttrium (Y)	[7440-65-5]	50	mg/L		
	Rhenium (Re)	[7440-15-5]	6,25	mg/L		
	Iridium (Ir)	[7439-88-5]	6,25	mg/L		
	Lithium (Li)	[7439-93-2]	35	mg/L		