

Diesel - Density - 833,0 kg/m3

Art. ID ROF-DE-GO-597
Unit 250 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1202 III E2

Description

The certified value of this CRM can be used with the following test methods: ISO 12185, ASTM D 4052, IP 365, ASTM D 1298, IP 160, ISO 3675

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Density		833	kg/m3		