

Custom Solution - 10 components 500 mg/L each of in HNO3 10%

Art. ID LM24-CUS-72518-2X100ML
Unit 2 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	500	mg/L		
	Calcium (Ca)	[7440-70-2]	500	mg/L		
	Copper (Cu)	[7440-50-8]	500	mg/L		
	Iron (Fe)	[7439-89-6]	500	mg/L		
	Magnesium (Mg)	[7439-95-4]	500	mg/L		
	Sodium (Na)	[7440-23-5]	500	mg/L		
	Nickel (Ni)	[7440-02-0]	500	mg/L		
	Phosphorus (P)	[7723-14-0]	500	mg/L		
	Vanadium (V)	[7440-62-2]	500	mg/L		
	Zinc (Zn)	[7440-66-6]	500	mg/L		