

Alloy Steel, Welding Wire, D6-AC, K24728, Ø 38 mm x 19 mm

Art. ID IARM-FED6AC-22-G
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified Value	Aluminium (Al)	[7429-90-5]	0,035	%		
Certified Value	Arsenic (As)	[7440-38-2]	0,0031	%		
Certified Value	Boron (B)	[7440-42-8]	0,00023	%		
Certified Value	Carbon (C)	[7440-44-0]	0,471	%		
Certified Value	Chromium (Cr)	[7440-47-3]	1,02	%		
Certified Value	Cobalt (Co)	[7440-48-4]	0,005	%		
Certified Value	Copper (Cu)	[7440-50-8]	0,036	%		
Certified Value	Manganese (Mn)	[7439-96-5]	0,82	%		
Certified Value	Molybdenum (Mo)	[7439-98-7]	0,98	%		
Certified Value	Nickel (Ni)	[7440-02-0]	0,539	%		
Certified Value	Niobium (Nb)	[7440-03-1]	0,0037	%		
Certified Value	Nitrogen (N)	[7727-37-9]	0,0038	%		
Certified Value	Oxygen (O)	[7782-44-7]	0,001	%		
Certified Value	Phosphorus (P)	[7723-14-0]	0,0068	%		
Certified Value	Silicon (Si)	[7440-21-3]	0,26	%		
Certified Value	Tin (Sn)	[7440-31-5]	0,0019	%		
Certified Value	Titanium (Ti)	[7440-32-6]	0,0021	%		
Certified Value	Tungsten (W)	[7440-33-7]	0,0026	%		
Certified Value	Vanadium (V)	[7440-62-2]	0,108	%		
Indicative Value	Lead (Pb)	[7439-92-1]	~20	ppm		
Indicative Value	Sulfur (S)	[7704-34-9]	~19	ppm		