

Trace Elements in Water (lot 0425)

Art. ID NW-TM-25.7
Unit 250 mL
Deliverydetails No Dangerous Good /not restricted

Description

Expiry date (lot 0820): 12/2023

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	29,4 ± 4,0	µg/L		
	Antimony (Sb)	[7440-36-0]	23,5 ± 2,4	µg/L		
	Arsenic (As)	[7440-38-2]	27,8 ± 2,3	µg/L		
	Barium (Ba)	[7440-39-3]	41,2 ± 3,4	µg/L		
	Beryllium (Be)	[7440-41-7]	26,8 ± 3,2	µg/L		
	Bismuth (Bi)	[7440-69-9]	18,8 ± 2,6	µg/L		
	Cadmium (Cd)	[7440-43-9]	23,3 ± 1,7	µg/L		
	Chromium (Cr)	[7440-47-3]	23,6 ± 1,5	µg/L		
	Cobalt (Co)	[7440-48-4]	27,0 ± 1,9	µg/L		
	Copper (Cu)	[7440-50-8]	25,2 ± 2,0	µg/L		
	Gallium (Ga)	[7440-55-3]	0,112 ± 0,025	µg/L		
	Iron (Fe)	[7439-89-6]	31,1 ± 4,0	µg/L		
	Lead (Pb)	[7439-92-1]	26,0 ± 2,9	µg/L		
	Lithium (Li)	[7439-93-2]	24,4 ± 2,6	µg/L		
	Manganese (Mn)	[7439-96-5]	25,1 ± 1,8	µg/L		
	Molybdenum (Mo)	[7439-98-7]	26,9 ± 2,2	µg/L		
	Nickel (Ni)	[7440-02-0]	15,6 ± 1,3	µg/L		
	Rubidium (Rb)	[7440-17-7]	18,7 ± 1,1	µg/L		
	Selenium (Se)	[7782-49-2]	29,1 ± 3,2	µg/L		
	Strontium (Sr)	[7440-24-6]	74,2 ± 5,2	µg/L		
	Thallium (Tl)	[7440-28-0]	31,7 ± 2,6	µg/L		
	Tin (Sn)	[7440-31-5]	23,7 ± 2,3	µg/L		
	Titanium (Ti)	[7440-32-6]	24,8 ± 2,2	µg/L		
	Uranium (U)	[7440-61-1]	26,3 ± 2,4	µg/L		
	Vanadium (V)	[7440-62-2]	27,3 ± 1,9	µg/L		
	Zinc (Zn)	[7440-66-6]	45,6 ± 5,6	µg/L		