

Custom Solution - 5 components 100 mg/L each of in HNO3 2%

Art. ID LM24-CUS-70886-1X100ML
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Gallium (Ga)	[7440-55-3]	100	mg/L		
	Yttrium (Y)	[7440-65-5]	100	mg/L		
	Indium (In)	[7440-74-6]	100	mg/L		
	Terbium (Tb)	[7440-27-9]	100	mg/L		
	Bismuth (Bi)	[7440-69-9]	100	mg/L		