

Human Serum (anti-tTG IgA and IgG)

Art. ID ERM-DA487/IFCC

Unit vial

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Anti-tissue transglutam inase (anti-tTG) IgA an tibodies	anti-tTG IgA antibodies		400 ± 18	IU/mL		
Anti-tissue transglutam inase (anti-tTG) IgG an tibodies	anti-tTG IgG antibodies		200 ± 25	IU/mL		