

Custom Solution - 9 components 100 mg/L each of in HNO3 5%

Art. ID LM24-CUS-73050
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Vanadium (V)	[7440-62-2]	100	mg/L		
	Aluminium (Al)	[7429-90-5]	100	mg/L		
	Barium (Ba)	[7440-39-3]	100	mg/L		
	Calcium (Ca)	[7440-70-2]	100	mg/L		
	Cerium (Ce)	[7440-45-1]	100	mg/L		
	Iron (Fe)	[7439-89-6]	100	mg/L		
	Molybdenum (Mo)	[7439-98-7]	100	mg/L		
	Titanium (Ti)	[7440-32-6]	100	mg/L		
	Magnesium (Mg)	[7439-95-4]	100	mg/L		