

Custom Solution - 3 components 100 mg/L each of in HNO3 5%, HF 0.1%

Art. ID LM24-CUS-73051
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) II E2

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tin (Sn)	[7440-31-5]	100	mg/L		
	Iron (Fe)	[7439-89-6]	100	mg/L		
	Molybdenum (Mo)	[7439-98-7]	100	mg/L		