

Lithium bis(trifluoromethanesulfonyl)imide (TFSI) (unlabeled) 100 µg/mL in Methanol

Art. ID CIL-ULM-11742-1.2

Unit 1,2 mL

Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lithium bis(trifluoromethanesulfonyl)imide (TFSI)	[90076-65-6]	100	µg/mL		