

cis-Permethrin (phenoxy-13C6, 99%) 50 µg/mL in Acetonitrile

Art. ID CIL-CLM-7322-A-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	cis-Permethrin (phenoxy -13C6,99%)	[61949-76-6]/un labeled	50	µg/mL		