

Custom Solution - 18 components in HNO₃ 5%, HCl 5%

Art. ID LM24-CUS-73060

Unit 1 x 250 mL bottle

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	20	mg/L		
	Calcium (Ca)	[7440-70-2]	10	mg/L		
	Cerium (Ce)	[7440-45-1]	100	mg/L		
	Copper (Cu)	[7440-50-8]	100	mg/L		
	Europium (Eu)	[7440-53-1]	100	mg/L		
	Iron (Fe)	[7439-89-6]	100	mg/L		
	Indium (In)	[7440-74-6]	100	mg/L		
	Potassium (K)	[7440-09-7]	100	mg/L		
	Lithium (Li)	[7439-93-2]	20	mg/L		
	Manganese (Mn)	[7439-96-5]	50	mg/L		
	Nickel (Ni)	[7440-02-0]	100	mg/L		
	Phosphorus (P)	[7723-14-0]	100	mg/L		
	Sulfur (S)	[7704-34-9]	500	mg/L		
	Scandium (Sc)	[7440-20-2]	50	mg/L		
	Strontium (Sr)	[7440-24-6]	20	mg/L		
	Vanadium (V)	[7440-62-2]	100	mg/L		
	Yttrium (Y)	[7440-65-5]	100	mg/L		
	Sodium (Na)	[7440-23-5]	50	mg/L		