

Custom Solution - 4 components in HNO3 5%, HCl 5%, HF 0.2%

Art. ID LM24-CUS-73061

Unit 1 x 250 mL bottle

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Molybdenum (Mo)	[7439-98-7]	50	mg/L		
	Titanium (Ti)	[7440-32-6]	100	mg/L		
	Zirconium (Zr)	[7440-67-7]	100	mg/L		
	Silicon (Si)	[7440-21-3]	100	mg/L		