

Custom Solution - 10 components in HNO₃ 2%

Art. ID LM24-CUS-73064
Unit 5 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	15	mg/L		
	Cobalt (Co)	[7440-48-4]	15	mg/L		
	Copper (Cu)	[7440-50-8]	100	mg/L		
	Iron (Fe)	[7439-89-6]	15	mg/L		
	Nickel (Ni)	[7440-02-0]	100	mg/L		
	Lead (Pb)	[7439-92-1]	100	mg/L		
	Tin (Sn)	[7440-31-5]	100	mg/L		
	Titanium (Ti)	[7440-32-6]	2	mg/L		
	Thallium (Tl)	[7440-28-0]	2	mg/L		
	Zinc (Zn)	[7440-66-6]	2	mg/L		