

Inorganics in Geoduck Clam Tissue (*Panopea generosa*)

Art. ID	NIST-2983
Unit	12 g
Deliverydetails	Dry ice shipment /not restricted

Description

The certified values delivered by this Standard Reference Material (SRM) are intended for validating methods for determining total arsenic (As) and inorganic arsenic (iAs) in geoduck clam tissue and similar materials. A unit of SRM 2983 consists of one jar sealed inside an aluminized pouch. The jar contains approximately 12 g of cryogenically homogenized frozen tissue material. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified Mass Fraction	Total arsenic		3,677 ± 0,063	mg/kg		
Value						
Certified Mass Fraction	Inorganic arsenic (iAs)		0,193 ± 0,011	mg/kg		
Value						