

Custom Solution - 24 components in HNO3 5%

Art. ID LM24-CUS-73077
Unit 1 x 500 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 8 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	0,2	mg/L		
	Calcium (Ca)	[7440-70-2]	0,2	mg/L		
	Magnesium (Mg)	[7439-95-4]	0,2	mg/L		
	Manganese (Mn)	[7439-96-5]	0,2	mg/L		
	Strontium (Sr)	[7440-24-6]	0,2	mg/L		
	Titanium (Ti)	[7440-32-6]	0,2	mg/L		
	Aluminium (Al)	[7429-90-5]	0,5	mg/L		
	Arsenic (As)	[7440-38-2]	0,5	mg/L		
	Cadmium (Cd)	[7440-43-9]	0,5	mg/L		
	Chromium (Cr)	[7440-47-3]	0,5	mg/L		
	Copper (Cu)	[7440-50-8]	0,5	mg/L		
	Iron (Fe)	[7439-89-6]	0,5	mg/L		
	Lithium (Li)	[7439-93-2]	0,5	mg/L		
	Molybdenum (Mo)	[7439-98-7]	0,5	mg/L		
	Nickel (Ni)	[7440-02-0]	0,5	mg/L		
	Ruthenium (Ru)	[7440-18-8]	0,5	mg/L		
	Tin (Sn)	[7440-31-5]	0,5	mg/L		
	Zinc (Zn)	[7440-66-6]	0,5	mg/L		
	Cobalt (Co)	[7440-48-4]	1	mg/L		
	Lead (Pb)	[7439-92-1]	1	mg/L		
	Antimony (Sb)	[7440-36-0]	1	mg/L		
	Boron (B)	[7440-42-8]	5	mg/L		
	Bismuth (Bi)	[7440-69-9]	5	mg/L		
	Silicon (Si)	[7440-21-3]	5	mg/L		