

Custom Solution - 16 components 100 mg/L in HNO3 5%

Art. ID LM24-CUS-73083
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	100	mg/L		
	Calcium (Ca)	[7440-70-2]	100	mg/L		
	Chromium (Cr)	[7440-47-3]	100	mg/L		
	Cobalt (Co)	[7440-48-4]	100	mg/L		
	Copper (Cu)	[7440-50-8]	100	mg/L		
	Iron (Fe)	[7439-89-6]	100	mg/L		
	Magnesium (Mg)	[7439-95-4]	100	mg/L		
	Manganese (Mn)	[7439-96-5]	100	mg/L		
	Molybdenum (Mo)	[7439-98-7]	100	mg/L		
	Nickel (Ni)	[7440-02-0]	100	mg/L		
	Phosphorus (P)	[7723-14-0]	100	mg/L		
	Potassium (K)	[7440-09-7]	100	mg/L		
	Selenium (Se)	[7782-49-2]	100	mg/L		
	Sodium (Na)	[7440-23-5]	100	mg/L		
	Tin (Sn)	[7440-31-5]	100	mg/L		
	Zinc (Zn)	[7440-66-6]	100	mg/L		