

Carbon Steel, Rephosphorized and Resulfurized, AISI 12L14, G12144, Disc 38 mm

Art. ID IARM-Fe12L14-21-X
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,0013 ± 0,0005	%		
	Arsenic (As)	[7440-38-2]	0,0027 ± 0,0005	%		
	Bismuth (Bi)	[7440-69-9]	0,0024 ± 0,0005	%		
	Carbon (C)	[7440-44-0]	0,063 ± 0,002	%		
	Cobalt (Co)	[7440-48-4]	0,0035 ± 0,0006	%		
	Chromium (Cr)	[7440-47-3]	0,056 ± 0,002	%		
	Copper (Cu)	[7440-50-8]	0,0134 ± 0,0006	%		
	Manganese (Mn)	[7439-96-5]	1,02 ± 0,02	%		
	Molybdenum (Mo)	[7439-98-7]	0,0034 ± 0,0005	%		
	Nitrogen (N)	[7727-37-9]	0,006 ± 0,001	%		
	Nickel (Ni)	[7440-02-0]	0,0177 ± 0,0007	%		
	Oxygen (O)	[7782-44-7]	0,012 ± 0,002	%		
	Phosphorus (P)	[7723-14-0]	0,071 ± 0,004	%		
	Lead (Pb)	[7439-92-1]	0,29 ± 0,01	%		
	Sulfur (S)	[7704-34-9]	0,32 ± 0,01	%		
	Antimony (Sb)	[7440-36-0]	0,0033 ± 0,0006	%		
	Tin (Sn)	[7440-31-5]	0,0013 ± 0,0004	%		
	Vanadium (V)	[7440-62-2]	0,0032 ± 0,0004	%		
	Tungsten (W)	[7440-33-7]	0,0030 ± 0,0009	%		
	Zirconium (Zr)	[7440-67-7]	0,0009 ± 0,0003	%		