

### **Custom Solution - N(NO<sub>2</sub>-) 5 µg/L in Water**

Art. ID LM24-CUS-72171-1x500ML  
Unit 1 x 500 mL bottle  
Deliverydetails No Dangerous Good /not restricted

#### **Description**

Conservation 8 months (originally closed bottle).

| Text/Information | Analyte/Parameter                    | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--------------------------------------|------------|---------------------|------|--------|--------|
|                  | Nitrites as N (N(NO <sub>2</sub> -)) |            | 5                   | µg/L |        |        |