

Custom Solution - 7 components in HNO3 2%

Art. ID LM24-CUS-72198-1X250ML
Unit 1 x 250 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 8 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	50	µg/L		
	Arsenic (As)	[7440-38-2]	10	µg/L		
	Chromium (Cr)	[7440-47-3]	25	µg/L		
	Iron (Fe)	[7439-89-6]	50	µg/L		
	Manganese (Mn)	[7439-96-5]	20	µg/L		
	Nickel (Ni)	[7440-02-0]	20	µg/L		
	Lead (Pb)	[7439-92-1]	10	µg/L		