

PARC 10 PFAS Native Mixture in Methanol (with 4 molar equivalents NaOH)

Art. ID CIL-ES-5670

Unit 1,2 mL

Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sodium Trifluoroacetate	[2923-18-4]	2000	ng/mL		
	Pentafluoropropanoic acid (PFPrA)	[422-64-0]	2000	ng/mL		
	Perfluorobutyric Acid, Sodium Salt	[2218-54-4]	2000	ng/mL		
	Perfluoropentanoic acid (PFPeA), sodium salt	[2706-89-0]	2000	ng/mL		
	Perfluorohexanoic acid (PFHxA), sodium salt	[2923-26-4]	2000	ng/mL		
	Sodium trifluoromethane sulfonate	[2926-30-9]	2000	ng/mL		
	Sodium perfluoroethanesulfonate (PFETs)	[2923-21-9]	2000	ng/mL		
	Sodium perfluoro-1-propanesulfonate (PFPrS)	[359868-82-9]	2000	ng/mL		
	Potassium perfluoro-1-butanesulfonate	[29420-49-3]	2000	ng/mL		
	Perfluoropentanesulfonate (PFPeS), sodium salt	[630402-22-1]	2000	ng/mL		