

## **Custom Solution - Dibutyl sulfide as Sulphur [544-40-1] 300 mg/L in Iso-octane**

Art. ID LM24-CUS-73095  
Unit 1 x 200 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1262 II E2

### Description

Conservation 12 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Dibutyl sulfide as Sulp hur	Di-N-butyl Sulfide (as S)	[544-40-1]	300	mg/L		