

Custom Solution - 15 components in HNO3 5%

Art. ID LM24-CUS-73096

Unit 2 x 50 mL bottle

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	100	mg/L		
	Barium (Ba)	[7440-39-3]	2000	mg/L		
	Cadmium (Cd)	[7440-43-9]	2	mg/L		
	Cr6+		20	mg/L		
	Iron (Fe)	[7439-89-6]	200	mg/L		
	Manganese (Mn)	[7439-96-5]	200	mg/L		
	Nickel (Ni)	[7440-02-0]	200	mg/L		
	Lead (Pb)	[7439-92-1]	20	mg/L		
	Copper (Cu)	[7440-50-8]	100	mg/L		
	Tin (Sn)	[7440-31-5]	1000	mg/L		
	Zinc (Zn)	[7440-66-6]	50	mg/L		
	Boron (B)	[7440-42-8]	500	mg/L		
	Antimony (Sb)	[7440-36-0]	250	mg/L		
	Selenium (Se)	[7782-49-2]	25	mg/L		
	Strontium (Sr)	[7440-24-6]	1250	mg/L		