

**EasyMP™ fragments (10-100 µm) in Ethanol (20 vol.%) + Tween 20 (0.05 vol.%) solution - different polymer types and colors available - please contact us for details**

Art. ID                      MPS-019  
Unit                        30 mL  
Deliverydetails            ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

**Description**

available colors: Clear/Red/Blue/Violet/Fluorescent green (Violet/Green only for PE)

| Text/Information        | Analyte/Parameter                                                                                                                                                                                                                                                                                                                                                            | CAS number | Concentration/Value | Unit | Method | Source |
|-------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------------|------|--------|--------|
| available polymer types | PE (polyethylene) / PP (polypropylene) / PVC (polyvinyl chloride) / PET (polyethylene terephthalate) / PS (polystyrene) / PUR (polyurethane) / TPU (thermoplastic polyurethane) / PS (polystyrene) / PA6 (Polyamide 6) / PA12 (Polyamide 12) / PC (polycarbonate) / CA (cellulose acetate) / PLA (polylactic acid) / PTFE (polytetrafluoroethylene) / POM (polyoxymethylene) |            | <1                  | %    |        |        |