

EasyMP™ fragments (100-500 µm) in Ethanol (20 vol.%) + Tween 20 (0.05 vol.%) solution - different polymer types and colors available - please contact us for details

Art. ID MPS-021
Unit 30 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

Description

available colors: Clear/Black/Red/Blue

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
available polymer types	PE (polyethylene) / PP (polypropylene) / PVC (polyvinyl chloride) / PET (polyethylene terephthalate) / PUR (polyurethane) / TPU (thermoplastic polyurethane) / PS (polystyrene) / ABS (acrylonitrile butadiene styrene) / PA6 (Polyamide 6) / PA12 (Polyamide 12) / PMMA (Polymethyl methacrylate) / PC (polycarbonate) / CA (cellulose acetate) / PLA (polylactic acid) / POM (polyoxymethylene)		<1	%		