

EasyMP™ beads (10-30 µm) in Ethanol (20 vol.%) + Tween 20 (0.05 vol.%) solution - different polymer types and colors available - please contact us for details

Art. ID MPS-025
Unit 30 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

Description

available colors: Clear/Red/Fluorescent green (Color only for PE)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
available polymer types	PE (polyethylene) / PS (polystyrene) / PMMA (P olymethyl methacrylate)		<1	%		