

EasyMP™ beads (40-60 µm) in Ethanol (20 vol.%) + Tween 20 (0.05 vol.%) solution - different polymer types and colors available - please contact us for details

Art. ID MPS-028
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG)

Description

available colors: Clear/Red/Blue/Fluorescent green (Color only for PE)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
available polymer types	PE (polyethylene) / PMMA (Polymethyl methacrylate)		<1	%		