

Custom Solution - 18 components 10 mg/L each of in HNO3 2%

Art. ID LM24-CUS-73114
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	10	mg/L		
	Boron (B)	[7440-42-8]	10	mg/L		
	Barium (Ba)	[7440-39-3]	10	mg/L		
	Calcium (Ca)	[7440-70-2]	10	mg/L		
	Chromium (Cr)	[7440-47-3]	10	mg/L		
	Caesium (Cs)	[7440-46-2]	10	mg/L		
	Copper (Cu)	[7440-50-8]	10	mg/L		
	Iron (Fe)	[7439-89-6]	10	mg/L		
	Potassium (K)	[7440-09-7]	10	mg/L		
	Lithium (Li)	[7439-93-2]	10	mg/L		
	Magnesium (Mg)	[7439-95-4]	10	mg/L		
	Sodium (Na)	[7440-23-5]	10	mg/L		
	Nickel (Ni)	[7440-02-0]	10	mg/L		
	Lead (Pb)	[7439-92-1]	10	mg/L		
	Rubidium (Rb)	[7440-17-7]	10	mg/L		
	Silicon (Si)	[7440-21-3]	10	mg/L		
	Ytterbium (Yb)	[7440-64-4]	10	mg/L		
	Zinc (Zn)	[7440-66-6]	10	mg/L		