

**CHEMPLEX, PreAssembled Sample Cup, Stackable Analyzer TrimLess® XRF Sample  
Cups, Single Open End, 1.69" (42.9mm) O.D. X 1.84" (46.7mm) Max O.D. X 0.73" (46.7mm)  
High with 12.0um Polypropylene - Designed to fit Horiba 6100/6800/60A Series**

|                 |                   |
|-----------------|-------------------|
| Art. ID         | CX-10P-8723       |
| Unit            | 22 sets/tube      |
| Deliverydetails | No Dangerous Good |