

**CHEMPLEX, PreAssembled Sample Cup, Stackable SpectroCup® with TrimLess® Sleeve,  
Internal Overflow Reservoir with Vented Cap and integrated knob, Double Open Ends,  
1.76" (44.7mm) O.D. X 1.85" (46.9mm) Max O.D. X 1.57" (39.9mm) High with 6.0um Mylar®**

|                 |                   |
|-----------------|-------------------|
| Art. ID         | CX-10P-P216       |
| Unit            | 10 sets/tube      |
| Deliverydetails | No Dangerous Good |