

Custom Solution - 6 components 10 mg/L each of in HNO3 5%/HF 0.2%

Art. ID LM24-CUS-73140
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Germanium (Ge)	[7440-56-4]	10	mg/L		
	Niobium (Nb)	[7440-03-1]	10	mg/L		
	Tantalum (Ta)	[7440-25-7]	10	mg/L		
	Titanium (Ti)	[7440-32-6]	10	mg/L		
	Tungsten (W)	[7440-33-7]	10	mg/L		
	Zirconium (Zr)	[7440-67-7]	10	mg/L		