

## Thermally refined Lead - chips

Art. ID                      BCR-287B  
Unit                         160 g  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,00152	%		
	Arsenic (As)	[7440-38-2]	<0,0000003	%		
	Bismuth (Bi)	[7440-69-9]	0,00673	%		
	Cadmium (Cd)	[7440-43-9]	0,000036	%		
	Copper (Cu)	[7440-50-8]	0,000098	%		
	Nickel (Ni)	[7440-02-0]	0,0000024	%		
	Antimony (Sb)	[7440-36-0]	0,000004	%		
	Selenium (Se)	[7782-49-2]	<0,000005	%		
	Tin (Sn)	[7440-31-5]	<0,000005	%		
	Tellurium (Te)	[13494-80-9]	<0,00002	%		
	Thallium (Tl)	[7440-28-0]	0,000073	%		
	Zinc (Zn)	[7440-66-6]	<0,00001	%		