

Nimonic 75 - Creep rate

Art. ID BCR-425
Unit set(3)
Deliverydetails No Dangerous Good /not restricted

Description

150 mm long 14 mm diameter bars of Nimonic 75. Testing conditions: T = 600 °C, s = 160 Mpa

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Creep rate at 400 h		72 x 10 ⁻⁶ /h			
	Time to 2 % strain		278 h			
	Time to 4 % strain		557 h			