

Estuarine water - Trace elements

Art. ID	BCR-505
Unit	1 L
Deliverydetails	No Dangerous Good /not restricted

Description

The material has been filtered at 0.45 µm and acidified with HNO₃ to around pH 1.5

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cadmium (Cd)	[7440-43-9]	0,8	nmol/kg		
	Copper (Cu)	[7440-50-8]	29,4	nmol/kg		
	Nickel (Ni)	[7440-02-0]	24,1	nmol/kg		
	Zinc (Zn)	[7440-66-6]	172	nmol/kg		