

Welding dust loaded on filter

Art. ID BCR-545
Unit ea
Deliverydetails No Dangerous Good /not restricted

Description

Welding dust loaded on a filter.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--------------------|--------------|---------------------|------|--------|--------|
| 40,2 | Chromium-VI (Cr+6) | [18540-29-9] | 40,2 | g/kg | | |
| 39,5 | total leachable Cr | | 39,5 | g/kg | | |