

Bovine eye (reconstituted) - Clenbuterol (blank)

Art. ID	BCR-673
Unit	vial
Deliverydetails	No Dangerous Good /not restricted

Description

The material consists of the equivalent of 2.00 ± 0.01 g (0.5 %) wet bovine eye material in sealed glass vials. It is a greyish powder obtained after lyophilisation and has to be reconstituted before use.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Mass fraction ($\mu\text{g}/\text{kg}$ wet eye material)	Clenbuterol free base		<0,50	$\mu\text{g}/\text{kg}$		