

**Bovine eye (reconstituted) - Clenbuterol (positive)**

Art. ID	BCR-674
Unit	vial
Deliverydetails	No Dangerous Good /not restricted

**Description**

The material consists of the equivalent of  $2.00 \pm 0.01$  g (0.5 %) wet bovine eye material in sealed glass vials. It is a greyish powder obtained after lyophilisation and has to be reconstituted before use.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Mass fraction ( $\mu\text{g}/\text{kg}$ wet eye material)	Clenbuterol free base		9	$\mu\text{g}/\text{kg}$		