

Linde type A zeolite

| | |
|-----------------|-----------------------------------|
| Art. ID | BCR-705 |
| Unit | 10 g |
| Deliverydetails | No Dangerous Good /not restricted |

Description

Adsorption of argon at the temperature of liquid argon (87 K) on a microporous material (Linde type A zeolite) according to DIN 66135-4.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|------------------------|------------|---------------------|--------------------|--------|--------|
| | Micropore volume | | 0.181 | cm ³ /g | | |
| | Median micropore width | | 0,592 | nm | | |