

## **Polychlorinated biphenyls in 2,2,4-Trimethylpentane (Isooctane)**

Art. ID                      ERM-AC823  
Unit                         1,2 mL  
Deliverydetails            ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

### Description

Indicative values for PCB 31, 77, 100, 149, 163, 170, 187 and 194

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,4,4'-Trichlorobiphenyl (PCB 28)	[7012-37-5]	703	µg/kg		
	2,2',5,5'-Tetrachlorobiphenyl (PCB 52)	[35693-99-3]	706	µg/kg		
	2,2',4,5,5'-Pentachlorobiphenyl (PCB 101)	[37680-73-2]	696	µg/kg		
	2,3',4,4',5-Pentachlorobiphenyl (PCB 118)	[31508-00-6]	712	µg/kg		
	2,2',3,4,4',5'-Hexachlorobiphenyl (PCB 138)	[35065-28-2]	678	µg/kg		
	2,2',4,4',5,5'-Hexachlorobiphenyl (PCB 153)	[35065-27-1]	702	µg/kg		
	2,2',3,4,4',5,5'-Heptachlorobiphenyl (PCB 180)	[35065-29-3]	700	µg/kg		