

BCR-ABL pDNA Calibrant (set of 6 calibrants to monitor chronic myelogenous leukemia)

Art. ID ERM-AD623
Unit set(6)
Deliverydetails Dry ice shipment

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ERM-AD623a: Copy number concentration of the p			1080000	cp/ μ L		
lasmid						
ERM-AD623b: Copy number concentration of the p			108000	cp/ μ L		
lasmid						
ERM-AD623c: Copy number concentration of the p			10300	cp/ μ L		
lasmid						
ERM-AD623d: Copy number concentration of the p			1020	cp/ μ L		
lasmid						
ERM-AD623e: Copy number concentration of the p			104	cp/ μ L		
lasmid						
ERM-AD623f: Copy number concentration of the p			10	cp/ μ L		
lasmid						