

**BCR-ABL pDNA Calibrant (set of 6 calibrants to monitor chronic myelogeneous leukemia)**

Art. ID                      ERM-AD623  
Unit                         set(6)  
Deliverydetails            Dry ice shipment

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ERM-AD623a: Copy number concentration of the p lasmid			1080000	cp/μL		
ERM-AD623b: Copy number concentration of the p lasmid			108000	cp/μL		
ERM-AD623c: Copy number concentration of the p lasmid			10300	cp/μL		
ERM-AD623d: Copy number concentration of the p lasmid			1020	cp/μL		
ERM-AD623e: Copy number concentration of the p lasmid			104	cp/μL		
ERM-AD623f: Copy number concentration of the p lasmid			10	cp/μL		