

**BCR-ABL pDNA Calibrant (set of 6 calibrants to monitor chronic myelogeneous leukemia)**

Art. ID                      ERM-AD623  
Unit                         set(6)  
Deliverydetails            Dry ice shipment

| Text/Information  | Analyte/Parameter | CAS number | Concentration/Value | Unit  | Method | Source |
|---|-------------------|------------|---------------------|-------|--------|--------|
| ERM-AD623a: Copy number<br>concentration of the p<br>lasmid |                   |            | 1080000             | cp/μL |        |        |
| ERM-AD623b: Copy number<br>concentration of the p<br>lasmid |                   |            | 108000              | cp/μL |        |        |
| ERM-AD623c: Copy number<br>concentration of the p<br>lasmid |                   |            | 10300               | cp/μL |        |        |
| ERM-AD623d: Copy number<br>concentration of the p<br>lasmid |                   |            | 1020                | cp/μL |        |        |
| ERM-AD623e: Copy number<br>concentration of the p<br>lasmid |                   |            | 104                 | cp/μL |        |        |
| ERM-AD623f: Copy number<br>concentration of the p<br>lasmid |                   |            | 10                  | cp/μL |        |        |