

Stacked Genetically modified MON 863 x MON 810 Maize (blank)

Art. ID ERM-BF417A
Unit vial
Deliverydetails No Dangerous Good /not restricted

Description

No contamination was detected in the non-GM material when applying event-specific MON 863 realtime PCR with a detection limit of 0.8 g/kg and no contamination was detected when applying eventspecific MON 810 real-time PCR with a detection limit of 1.1g /kg. With a confidence level of 95 % the MON 863 x MON 810 maize mass fraction is below the certified value.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-----------------------------|------------|---------------------|------|--------|--------|
| | MON 863 x MON 810 maiz e | | <1,0 | g/kg | | |