

Genetically modified H7-1 Sugar Beet (level 1)

Art. ID	ERM-BF419B
Unit	vial
Deliverydetails	No Dangerous Good /not restricted

Description

ERM-BF419b is available in glass bottles containing approximately 1 g of sugar beet powder closed under argon atmosphere. This reference material has been produced from genetically modified sugar beet roots delivered by KWS (Einbeck, DE). According to Commission Regulation (EC) No 65/2004 the event H7-1 sugar beet corresponds to the unique identifier KM-???H71-4. The certified H7-1 mass fraction is based on the mass fraction of modified sugar beet roots used for the processing of the CRM. All sugar beet roots used in the processing were individually tested by ELISA and event-specific rt-PCR.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	H7-1 sugar beet		1000	g/kg		