

Chloroform-d,"100%", 99.96 atom % D, contains 0.5 wt. % silver wire as stabilizer (ALDRICH)

Art. ID SAF-431915-50ML
Unit 50 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1888 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| | Chloroform-d1 | [865-49-6] | | | | |