

Cortisol in human serum (unspiked)

Art. ID ERM-DA192
Unit ampoule
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Amount of cortisol in the ampoule	[50-23-7]	0,1315 ± 0,0025	µg		
	Amount of cortisol in the ampoule	[50-23-7]	0,363 ± 0,007	nmol		
	Concentration of cortisol in the reconstituted serum	[50-23-7]	98,8 ± 2,0	µg/L		
	Concentration of cortisol in the reconstituted serum	[50-23-7]	273 ± 6	nmol/L		