

Cortisol in human serum (spiked)

Art. ID ERM-DA193
Unit ampoule
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Amount of cortisol in the ampoule	[50-23-7]	0,370 ± 0,003	µg		
	Amount of cortisol in the ampoule	[50-23-7]	1,022 ± 0,008	nmol		
	Concentration of cortisol in the reconstituted serum	[50-23-7]	277 ± 5	µg/L		
	Concentration of cortisol in the reconstituted serum	[50-23-7]	763 ± 14	nmol/L		