

n-Nonane - Flash point

Art. ID	ERM-FC032
Unit	100 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1920 III E1

Description

The certified value was determined by the Abel closed cup method described in the Institute of Petroleum Standard IP170/95 and also published as British Standard BS2000:Part 170: 1995. The certified value is corrected to standard barometric pressure at 0°C.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Non-equilibrium flash point		32,5 ± 0,5	°C		